


Disposable Auto Injector



Device Type

Disposable Auto Injector

Add technology-platform Container

1 ml to 1.5 ml standard cartridge

Dosage Type

Fixed Dose

Route of Administration

Subcutaneous

Needle Insertion

Manual/Automatic


Activation Method

Sleeve Press

Injection Feedback

Click at the start of injection, Dual visual feedback & Viewing window

Handling Instructions



Key Features

- ✓ Simple 3 step operation
- ✓ Sleeve activation
- ✓ Auto/manual needle insertion
- ✓ Initiation audible click and dual visual feedback
- ✓ One-handed operation
- ✓ Permanently hidden needle
- ✓ Passive locking needle cover
- ✓ Preconfigured or customize client specifications