





<b>Device Type</b>	1ml Glass PFS NSD
<b>Add technology-platform Container</b>	1ml Glass PFS
<b>Dosage Type</b>	Fixed Dose
<b>Route of Administration</b>	Subcutaneous
<b>Compatibility</b>	Compatible with marketed prefilled glass syringes
<b>Dose monitoring</b>	Wide window availability for easy monitoring
<b>Safety</b>	Fail safe & Robust needle coverage

## Handling Instructions:




*Remove The Cap*

2



*Inject the medicine dose in the human body  
by pressing thumb pad*

+



*Retract the device  
&  
actuate the NDL Cover*

### Key Features

- ✓ Simple 3 step operation
- ✓ Allows smooth operation
- ✓ self locking of cartridge after needle cover actuation injection
- ✓ Robust and ergonomic design
- ✓ Easy syringe assembly