



**Device Type**

Micro Pre-Filled Syringe - Disposable

**Container**

0.2 to 0.5ml

**Dosage Type**

Fixed Dose

**Route of Administration**

Hypodermic

**Needle Insertion**

Manual

**Activation Method**

Thumb pad

**Injection Feedback**

Visual monitoring window


## Handling Instructions:

1



*Remove The Cap*

2



*Penetrate the needle*

3



*Thumb pad pressing*

+



*Discard the device in sharps container*

### Key Features:

- ✓ Simple one step operation
- ✓ Plunger rod self locking after medicine dose injection
- ✓ Visual monitoring throughout the medicine dose injection

**List of drugs that can be used with micro prefilled syringe:**

- Clexane (Enoxaparin)
- Liraglutide Depot
- Humira
- Fraxiparine
- Fraxiparine Forte
- Aranesp
- Methotrexate

## Advantages:

- Ease of administration,
- More convenient for healthcare professional and end user,
- Easier for home use, easier in emerging situations,
- Reduction of medication errors,
- Misidentification, and better dose accuracy,
- Increase assurance of sterility,
- Better use of controlled drugs such as narcotics,
- Lower injection cost by less preparation,
- Fever material,
- Easy storage and disposal,
- Elimination of vials overfill for products transferred to syringes for direct diluents,
- Removal of preservatives, i.e. thiomersal from vaccine formulation, product differentiation